

Anti-CRK antibody [3G11E8] (STJ97969)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/FC/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES



Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:200-1:1000 IF 1:200-1:1000 FC 1:200-1:400 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG2b
Molecular Weight	Calculated: 34kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	1398
Gene Symbol	CRK
UniProt ID	CRK_HUMAN
Specificity	Crk II Monoclonal Antibody detects endogenous levels of Crk II protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Flow cytometric analysis of HeLa cells using Crk II monoclonal antibody (blue) and negative control (red).



Immunohistochemistry analysis of paraffin-embedded intima cancer tissues with DAB staining using Crk II monoclonal antibody.



Immunofluorescence analysis of 3T3-L1 cells using Crk II monoclonal antibody (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Western blot analysis using Crk II monoclonal antibody against HEK293 (1) and CRK-hlgGfc transfected HEK293 (2) cell lysate.