

Anti-EPHA8 antibody (70-150 aa) [9A12D8] (STJ98030)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 111kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2046
Gene Symbol	EPHA8
UniProt ID	EPHA8_HUMAN
Immunogen Region	70-150 aa
Specificity	EphA8 Monoclonal Antibody detects endogenous levels of EphA8 protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis using EphA8 monoclonal antibody against truncated Trx-EphA8 recombinant protein (1) and truncated MBP-EphA8 (aa70-150) recombinant protein (2).