

Anti-NR3C1 antibody [6E6] (STJ98115)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/FC/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES





Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:200-1:1000 IF 1:200-1:1000 FC 1:200-1:400 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 86kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2908
Gene Symbol	NR3C1
UniProt ID	GCR_HUMAN
Specificity	GR Monoclonal Antibody detects endogenous levels of GR protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Flow cytometric analysis of K562 cells using GR monoclonal antibody (green) and negative control (purple).

Immunohistochemistry analysis of paraffin-embedded lung cancer tissues with DAB staining using GR monoclonal antibody.

Immunofluorescence analysis of PC-2 cells using GR monoclonal antibody (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin



Western blot analysis using GR monoclonal antibody against HeLa cell lysate.