

Anti-HFE antibody [3F1] (STJ98132)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 40kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION



Gene ID	3077
Gene Symbol	HFE
UniProt ID	HFE_HUMAN
Specificity	HFE Monoclonal Antibody detects endogenous levels of HFE protein.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Immunofluorescence analysis of HepG2 cells using HFE monoclonal antibody (green). Blue: DRAQ5 fluorescent DNA dye.



Western blot analysis using HFE monoclonal antibody against HEK293 (1) and HFE-hlgGc transfected HEK293 (2) cell lysate.