

Mouse Anti-Human IgG antibody {Unconjugated} (STJ98146)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Specificity	Human IgG Monoclonal Antibody detects Human IgG protein.
-------------	--

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis using Human IgG monoclonal antibody (lane 1, 2) and ALB monoclonal antibody (lane 3, 4) against human serum (lane 1, 3) and plasma (lane 2, 4).