

Anti-IKBKE antibody (1-257 aa) [6B4B5] (STJ98166)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 80kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9641
Gene Symbol	IKBKE
UniProt ID	IKKE_HUMAN
Immunogen Region	1-257 aa
Specificity	IKK Epsilon Monoclonal Antibody detects endogenous levels of IKK Epsilon protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis using IKK Epsilon monoclonal antibody against truncated IKK Epsilon recombinant protein (1) and full-length IKK Epsilon (aa1-716) -hlgGfc transfected COS7 cell lysate (2).