

Anti-KSHV ORF45 antibody [2D4A5] (STJ98209)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/ELISA
Host / Source	Mouse
Reactivity	Sheep

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ICC 1:200-1:1000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Specificity	KSHV ORF45 Monoclonal Antibody detects KSHV ORF45 protein.
-------------	--

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunocytochemistry analysis of TPA induced BCBL-1 cells (A) and uninduced BCBL-1 cells (B) using KSHV ORF45 monoclonal antibody with AEC staining.



Western blot analysis using KSHV ORF45 monoclonal antibody against KSHV ORF45 recombinant protein.