

## Anti-LCK antibody [8E5F9/4B2C7] (STJ98215)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                     |
| Conjugation         | Unconjugated                                                                                   |
| Purification        | Affinity purification                                                                          |
| Dilution Range      | WB 1:500-1:2000<br>ELISA 1:10000                                                               |
| Formulation         | Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.                           |
| Isotype             | IgG2a                                                                                          |
| Molecular Weight    | Calculated: 58kD                                                                               |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|             |                                                                   |
|-------------|-------------------------------------------------------------------|
| Gene ID     | <a href="#">3932</a>                                              |
| Gene Symbol | <a href="#">LCK</a>                                               |
| UniProt ID  | <a href="#">LCK_HUMAN</a>                                         |
| Specificity | Lck Monoclonal Antibody detects endogenous levels of Lck protein. |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis using Lck monoclonal antibody against MOLT-4 (1), CCRF-CEM (2), CCRF-HSB-2 (3) and Jurkat (4) cell lysate.