

Anti-NEUROG3 antibody [6G2] (STJ98273)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Antibody are purified by protein G affinity chromatography.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 50% Glycerol and 0.03% Sodium Azide.
Isotype	IgG1
Molecular Weight	Calculated: 23kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50674
Gene Symbol	NEUROG3
UniProt ID	NGN3_HUMAN
Specificity	Neurogenin 3 Monoclonal Antibody detects endogenous levels of Neurogenin 3 protein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis using Neurogenin 3 monoclonal antibody against full-length NGN3 (aa1-214) -hlgGfC transfected HEK293 cell lysate (1).