

Anti-NME1 antibody [4B2] (STJ98280)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/FC/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:200-1:1000 IF 1:200-1:1000 FC 1:200-1:400 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 17kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4830
Gene Symbol	NME1
UniProt ID	NDKA_HUMAN
Specificity	NM23-H1 Monoclonal Antibody detects endogenous levels of NM23-H1 protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Flow cytometric analysis of Jurkat cells using NM23-H1 monoclonal antibody (green) and negative control (purple).



Immunohistochemistry analysis of paraffin-embedded human Placenta tissues with AEC staining using NM23-H1 monoclonal antibody.



Immunofluorescence analysis of HeLa cells using NM23-H1 monoclonal antibody (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin. Blue: DRAQ5 fluorescent DNA dye.



Western blot analysis using NM23-H1 monoclonal antibody against NME1-hlgGf transfected HEK293 cell lysate.