

Anti-SLC22A1 antibody [2C5] (STJ98284)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/FC/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 FC 1:200-1:400 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 61kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6580
Gene Symbol	SLC22A1
UniProt ID	S22A1_HUMAN
Specificity	OCT1 Monoclonal Antibody detects endogenous levels of 41183 protein.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Flow cytometric analysis of Jurkat cells using OCT1 monoclonal antibody (green) and negative control (purple).



Western blot analysis using OCT1 monoclonal antibody against HEK293 (1) and SLC22A1-hlgGfC transfected HEK293 (2) cell lysate.