

Anti-PAWR antibody (1-330 aa) [3G9H7/4H12E9] (STJ98306)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 37kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5074
Gene Symbol	PAWR
UniProt ID	PAWR_HUMAN
Immunogen Region	1-330 aa
Specificity	PAR4 Monoclonal Antibody detects endogenous levels of PAR4 protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraffin-embedded human brain (left) and breast carcinoma (right), showing cytoplasmic and membrane localization with DAB staining using PAR4 monoclonal antibody.



Western blot analysis using PAR4 monoclonal antibody against full-length Trx-PAR4 recombinant protein (1) and HeLa cell lysate (2).