

Anti-ATP2C1 antibody [4G12] (STJ98326)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Mouse
Reactivity	Human/Monkey

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 101kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	27032
Gene Symbol	ATP2C1
UniProt ID	AT2C1_HUMAN
Specificity	PMR1 Monoclonal Antibody detects endogenous levels of PMR1 protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraffin-embedded human ovarian cancer (left) and breast cancer (right) tissues with DAB staining using PMR1 monoclonal antibody.



Western blot analysis using PMR1 monoclonal antibody against A431 (1), HeLa (2) and HEK293 (3) cell lysate.