

Anti-RAB10 antibody [4E2] (STJ98346)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ELISA
Host / Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 23kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10890
Gene Symbol	RAB10
UniProt ID	RAB10_HUMAN
Specificity	Rab 10 Monoclonal Antibody detects endogenous levels of Rab 10 protein.

ADDITIONAL INFORMATION

Note: STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of LOVO cells using Rab 10 monoclonal antibody (green). Blue: DRAQ5 fluorescent DNA dye.



Western blot analysis using Rab 10 monoclonal antibody against HeLa (1) and NIH/3T3 (2) cell lysate.