

Anti-RET antibody (896-1063 aa) [8D10C9] (STJ98358)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 124kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5979
Gene Symbol	RET
UniProt ID	RET_HUMAN
Immunogen Region	896-1063 aa
Specificity	Ret Monoclonal Antibody detects endogenous levels of Ret protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis using Ret monoclonal antibody against truncated RET recombinant protein (1) and RET (aa658-1063) -hlgGfC transfected CHO-K1 cell lysate (2).