

Anti-RTN3 antibody [1E11] (STJ98370)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:200-1:1000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10313
Gene Symbol	RTN3
UniProt ID	RTN3_HUMAN
Specificity	Rtn-3 Monoclonal Antibody detects endogenous levels of Rtn-3 protein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Confocal immunofluorescence analysis of HeLa (A), A431 (B) and THP-1 (C) cells using Rtn-3 monoclonal antibody (green). Blue: DRAQ5 fluorescent DNA dye.