

Anti-ZBTB7B antibody [7C12] (STJ98421)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:200-1:1000 IF 1:200-1:1000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 58kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	51043
Gene Symbol	ZBTB7B
UniProt ID	ZBT7B_HUMAN
Specificity	TH-POK Monoclonal Antibody detects endogenous levels of TH-POK protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraffin-embedded human Breast tissues with AEC staining using TH-POK monoclonal antibody.



Immunofluorescence analysis of HeLa cells using TH-POK monoclonal antibody (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin



Western blot analysis using TH-POK monoclonal antibody against HEK293 (1, 2) cell lysate.