

Anti-TRIM5 antibody [3B11H2] (STJ98429)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 56kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	85363
Gene Symbol	TRIM5
UniProt ID	TRIM5_HUMAN
Specificity	TRIM5 Alpha Monoclonal Antibody detects endogenous levels of TRIM5 Alpha protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraffin-embedded human metastatic adenocarcinoma (A) and stomach adenocarcinoma (B), showing cytoplasmic localization with AEC staining (A) and DAB staining (B) using TRIM5 Alpha monoclonal antibody.



Western blot analysis using TRIM5 Alpha monoclonal antibody against human breast carcinoma tissue lysate.