

Anti-ARRB2 antibody (C-Term) (STJ98468)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Pig

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000 IF 1:100-1:500
Formulation	Liquid in 0.1M Tris-Glycine (pH7.4, 150 mM NaCl) , 0.2% Sodium Azide and 50% Glycerol.
Molecular Weight	Calculated: 46kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	409
Gene Symbol	ARRB2
UniProt ID	ARRB2_HUMAN
Immunogen Region	C-Term
Specificity	Arrestin-Beta-2 Monoclonal Antibody detects endogenous levels of Arrestin-Beta-2 protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of HeLa cells using Arrestin-Beta-2 monoclonal antibody.



Western blot analysis using Arrestin-Beta-2 monoclonal antibody against HepG2 cell lysate.