

## Anti-FOXO1 antibody (C-Term) (STJ98492)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IF              |
| Host / Source | Mouse              |
| Reactivity    | Human/Mouse        |

### PRODUCT PROPERTIES


|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                     |
| Concentration       | 1 mg/mL                                                                                        |
| Conjugation         | Unconjugated                                                                                   |
| Purification        | Affinity purification                                                                          |
| Dilution Range      | WB 1:1000-1:2000<br>IF 1:100-1:500                                                             |
| Formulation         | Liquid in 0.1M Tris-Glycine (pH7.4, 150 mM NaCl) , 0.2% Sodium Azide and 50% Glycerol.         |
| Molecular Weight    | Calculated: 70kD                                                                               |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                  |                                                                         |
|------------------|-------------------------------------------------------------------------|
| Gene ID          | <a href="#">2308</a>                                                    |
| Gene Symbol      | <a href="#">FOXO1</a>                                                   |
| UniProt ID       | <a href="#">FOXO1_HUMAN</a>                                             |
| Immunogen Region | C-Term                                                                  |
| Specificity      | FoxO1A Monoclonal Antibody detects endogenous levels of FoxO1A protein. |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis using FoxO1A monoclonal antibody against A431, HeLa cell lysate.