

## Anti-HDAC3 antibody (STJ98496)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Dog/Pig/Rabbit

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000
Formulation	Liquid in 0.1M Tris-Glycine (pH7.4, 150 mM NaCl) , 0.2% Sodium Azide and 50% Glycerol.
Molecular Weight	Calculated: 49kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">8841</a>
Gene Symbol	<a href="#">HDAC3</a>
UniProt ID	<a href="#">HDAC3_HUMAN</a>
Specificity	HDAC3 Monoclonal Antibody detects endogenous levels of HDAC3 protein.

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis using HDAC3 monoclonal antibody against THP-1, LNCAP, MCF7, PC12, 293, Jurkat, 3T3, HeLa cell lysate.