

Anti-RAB2A antibody (C-Term) (STJ98545)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Cow/Chicken/Dog/Pig/Zebrafish

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000 IF 1:100-1:500
Formulation	Liquid in 0.1M Tris-Glycine (pH7.4, 150 mM NaCl) , 0.2% Sodium Azide and 50% Glycerol.
Molecular Weight	Calculated: 24kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5862
Gene Symbol	RAB2A
UniProt ID	RAB2A_HUMAN
Immunogen Region	C-Term
Specificity	Rab 2A Monoclonal Antibody detects endogenous levels of Rab 2A protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of HeLa cells using Rab 2A monoclonal antibody.



Western blot analysis using Rab 2A monoclonal antibody against HeLa cell lysate.