

Anti-UHRF1 antibody (N-Term) (STJ98567)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES



Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000 IF 1:100-1:500
Formulation	Liquid in 0.1M Tris-Glycine (pH7.4, 150 mM NaCl) , 0.2% Sodium Azide and 50% Glycerol.
Molecular Weight	Calculated: 90kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29128
Gene Symbol	UHRF1
UniProt ID	UHRF1_HUMAN
Immunogen Region	N-Term
Specificity	UHRF1 Monoclonal Antibody detects endogenous levels of UHRF1 protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western Blot analysis using UHRF1 Monoclonal Antibody against mouse heart, 293T, HeLa cell lysate, HeLa nuclear extract. Immunofluorescence analysis of HeLa cells using UHRF1 Monoclonal Antibody.