

Anti-USP7 antibody [3E4-D10-E8] (STJ99076)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Molecular Weight	Observed: 135kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7874
Gene Symbol	USP7
UniProt ID	UBP7_HUMAN
Specificity	This antibody detects endogenous levels of HAUSP/USP7 and does not cross-react with related proteins.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot detection of HAUSP/USP7 in HeLa and K562 cell lysates using HAUSP/USP7 mouse mAb (1:200 diluted) , with super ECL. Predicted band size: 135kDa. Observed band size: 135kDa.