

Anti-GTF2H1 antibody [3E8-D3-H8] (STJ99084)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/ICC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:300 ICC 1:50
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Molecular Weight	Observed: 62kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2965
Gene Symbol	GTF2H1
UniProt ID	TF2H1_HUMAN
Specificity	This antibody detects endogenous levels of GTF2H1 and does not cross-react with related proteins.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using GTF2H1 mouse mAb (dilution 1:50).



Immunoprecipitation analysis of HeLa cell lysates using GTF2H1 mouse mAb.



Western blot detection of GTF2H1 in HeLa, HeLa NE and Jurkat cell lysates using GTF2H1 mouse mAb (1:300 diluted). Predicted band size: 62kDa. Observed band size: 62kDa.