

Anti-IKZF1 antibody (C-Term) [1A12-F2-D8] (STJ99100)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Molecular Weight	Observed: 58kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10320
Gene Symbol	IKZF1
UniProt ID	IKZF1_HUMAN
Immunogen Region	C-Term
Specificity	This antibody detects endogenous levels of Ikaros (C-terminus) and does not cross-react with related proteins.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunoprecipitation analysis of HeLa cell lysate using Ikaros (C-terminus) mouse mAb.



Western blot detection of Ikaros (C-terminus) in Ramos and Jurkat cell lysates using Ikaros (C-terminus) mouse mAb (1:1000 diluted). Predicted band size: 58 kDa. Observed band size: 58 kDa.