

Anti-TAB1 antibody (N-Term) [3D10-G10-B6] (STJ99124)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Molecular Weight	Observed: 55kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10454
Gene Symbol	TAB1
UniProt ID	TAB1_HUMAN
Immunogen Region	N-Term
Specificity	This antibody detects endogenous levels of TAB1 (N-terminus) and does not cross-react with related proteins.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot detection of TAB1 (N-terminus) in K562 and Hela lysates using TAB1 (N-terminus) mouse mAb (1:1000 diluted). Predicted band size: 55KDa. Observed band size: 55KDa.