

Anti-CHRD1 antibody [3H1-F6-A10] (STJ99125)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000 IHC-P 1:300 ICC 1:300
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Molecular Weight	Observed: 52kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	91851
Gene Symbol	CHRD1
UniProt ID	CRDL1_HUMAN
Specificity	This antibody detects endogenous levels of CHRD1 and does not cross-react with related proteins.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot detection of CHRD1 in human eye and HeLa cell lysates using CHRD1 mouse mAb (1:1000 diluted). Predicted STJ99125



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight). 2, STJ99125



Immunofluorescence (IF): Immunocytochemistry stain of HeLa using CHRD1 mouse mAb (1:300). STJ99125