

## Anti-LIN28A antibody [2C1-F9-A2] (STJ99177)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB                 |
| Host / Source | Mouse              |
| Reactivity    | Human/Mouse        |

### PRODUCT PROPERTIES


|                     |                                                                                                                    |
|---------------------|--------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                         |
| Concentration       | 1 mg/mL                                                                                                            |
| Conjugation         | Unconjugated                                                                                                       |
| Purification        | The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range      | WB 1:1000                                                                                                          |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                            |
| Isotype             | IgG2a                                                                                                              |
| Molecular Weight    | Observed: 26kD                                                                                                     |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                     |

### TARGET INFORMATION


|             |                                                                                                   |
|-------------|---------------------------------------------------------------------------------------------------|
| Gene ID     | <a href="#">79727</a>                                                                             |
| Gene Symbol | <a href="#">LIN28A</a>                                                                            |
| UniProt ID  | <a href="#">LN28A_HUMAN</a>                                                                       |
| Specificity | This antibody detects endogenous levels of LIN28A and does not cross-react with related proteins. |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot detection of LIN28A in LIN28B and LIN28A recombinant antigen fragments and using LIN28A mouse mAb (1:1000 diluted).



Western blot detection of LIN28A in NTERA and F9 cell lysates using LIN28A mouse mAb (1:1000 diluted). Predicted band size: 26KDa. Observed band size: 26KDa.