

Anti-EGF antibody [7C11-G10-H8] (STJ99272)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution Range	ELISA 1:10000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Molecular Weight	Observed:
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1950
Gene Symbol	EGF
UniProt ID	EGF_HUMAN
Specificity	This antibody detects Human EGF.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Standard Curve for Human EGF: Capture Antibody Mouse mAb 201165-5E9 (5E9-A10-A6) to Human EGF at 4 µg/ml and Detector Antibody Mouse mAb 201165-7C11 (7C11-G10-H8) to Human EGF at 1 µg/ml.