

## Anti-HTRA2 antibody [1H8-C12-F8] (STJ99318)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Mouse/Hamster

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Molecular Weight	Observed: 36kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">27429</a>
Gene Symbol	<a href="#">HTRA2</a>
UniProt ID	<a href="#">HTRA2_HUMAN</a>
Specificity	This antibody detects endogenous levels of HtrA2/Omi and does not cross-react with related proteins.

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot detection of HtrA2/Omi in CHO-K1, MEF and 3T3 cell lysates using HtrA2/Omi mouse mAb (dilution 1:1000). Predicted band size: 49kDa. Observed band size: 36kDa.