

## Rabbit Anti-Horse IgG Heavy & Light Chain antibody {FITC} (STJ99461)

### GENERAL INFORMATION

|               |                      |
|---------------|----------------------|
| Product Type  | Secondary antibodies |
| Applications  | WB/IHC/IF/ELISA      |
| Host / Source | Rabbit               |
| Reactivity    | Horse                |

### PRODUCT PROPERTIES

|                     |                                                                                                 |
|---------------------|-------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                      |
| Concentration       | 2 mg/mL                                                                                         |
| Conjugation         | FITC                                                                                            |
| Purification        | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| Dilution Range      | IF 1:250-1:2000<br>FCM 1:250-1:2000                                                             |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                         |
| Isotype             | IgG                                                                                             |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.  |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|