

Rabbit Anti-Guinea Pig IgG Heavy & Light Chain antibody {FITC} (STJ99465)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Rabbit
Reactivity	Guinea Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	FITC
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Dilution Range	IF 1:250-1:2000 FCM 1:250-1:2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---