

## Mouse Anti-Human Ig Kappa Chain antibody {HRP} (STJ99535)

### GENERAL INFORMATION

|               |                      |
|---------------|----------------------|
| Product Type  | Secondary antibodies |
| Applications  | WB/IHC/IF/ELISA      |
| Host / Source | Mouse                |
| Reactivity    | Human                |

### PRODUCT PROPERTIES

|                     |                                                                                                 |
|---------------------|-------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                      |
| Concentration       | 2 mg/mL                                                                                         |
| Conjugation         | HRP                                                                                             |
| Purification        | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| Dilution Range      | WB 1:1000-10000<br>ELISA 1:10000-20000<br>IHC 1:50-100<br>IF 1:50-1000                          |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                         |
| Isotype             | Ig                                                                                              |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.  |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|