

Rabbit Anti-Human IgG Gamma chain antibody {Unconjugated} (STJ99586)

GENERAL INFORMATION

| | |
|---------------|----------------------|
| Product Type | Secondary antibodies |
| Applications | WB/IHC/IF/ELISA |
| Host / Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Polyclonal |
| Concentration | 2 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| Dilution Range | WB 1:1000-10000 ELISA 1:10000-20000 IHC 1:50-100 IF 1:50-1000 |
| Formulation | Liquid in PBS containing 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|