

Anti-MHC 2 antibody [M5/114] {PE/Cyanine5} (STJA0000048)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	14961
Gene Symbol	H2-Ab1
UniProt ID	HB2A_MOUSE

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Mouse splenocytes were stained with Anti-Mouse I-A/I-E Monoclonal Antibody (PE/Cyanine5 Conjugated) (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.