

## Anti-MHC 2 antibody [M5/114] {PE} (STJA0000054)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">14961</a>
Gene Symbol	<a href="#">H2-Ab1</a>
UniProt ID	<a href="#">HB2A_MOUSE</a>

### ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Mouse splenocytes were stained with Anti-Mouse I-A/I-E Monoclonal Antibody (PE Conjugated) [Used at 0.02 µg/10<sup>6</sup> cells dilution] (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.