

Anti-MHC 2 antibody [M5/114] {PerCP/Cyanine5.5} (STJA0000059)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PerCP/Cyanine5.5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	14961
Gene Symbol	H2-Ab1
UniProt ID	HB2A_MOUSE

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



C57BL/6 murine splenocytes were stained with Anti-Mouse I-A/I-E Monoclonal Antibody (PerCP/Cyanine5.5 Conjugated)[Used at 0.05 µg/10⁶ cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.