

Anti-CD11c/Ilgax antibody [N418] {FITC} (STJA0000062)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16411
Gene Symbol	Ilgax
UniProt ID	ITAX_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine splenocytes were stained with Anti-CD11c antibody (STJA0000062) and Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (APC Conjugated) (Left). Splenocytes stained with Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (APC Conjugated) (Right) were used as control