

Anti-CD11c/Ilgax antibody [N418] {PerCP/Cyanine5.5} (STJA0000068)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16411
Gene Symbol	Ilgax
UniProt ID	ITAX_MOUSE

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



C57BL/6 murine splenocytes were stained with Anti-Mouse CD11c Monoclonal Antibody (PerCP-Cy5.5 Conjugated) and Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (FITC Conjugated) (left), Splenocytes stained with Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (FITC Conjugated) (right) were used as control.