

## Anti-CD11c/Ilgax antibody [N418] {PE} (STJA0000072)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Armenian Hamster
Reactivity	Mouse

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">16411</a>
Gene Symbol	<a href="#">Ilgax</a>
UniProt ID	<a href="#">ITAX_MOUSE</a>

### ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine splenocytes were stained with Anti-CD11c antibody (STJA0000072) [Used at 0.02 µg/10<sup>6</sup> cells dilution] and Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (APC Conjugated). Splenocytes stained with Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (APC Conjugated) (Right) were used as control