

Anti-Cd80 antibody [16-10A1] {PE/Cyanine5} (STJA0000082)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12519
Gene Symbol	Cd80
UniProt ID	CD80_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Mouse abdominal macrophages elicited by starch broth were stained with Anti-CD80 antibody (STJA0000082) (filled gray histogram). Unstained macrophages (blank black histogram) were used as control.