

Anti-CD83 antibody [HB15e] {APC} (STJA0000093)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	APC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	9308
Gene Symbol	CD83
UniProt ID	CD83_HUMAN

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



PMA and ionomycin-stimulated (4h) Jurkat cells were stained with Anti-CD83 antibody (STJA0000093) (filled gray histogram) or Mouse IgG1, κ Isotype Control (APC Conjugated) (empty black histogram).