

Anti-Cd86 antibody [GL-1] {FITC} (STJA0000104)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12524
Gene Symbol	Cd86
UniProt ID	CD86_MOUSE

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



LPS-stimulated (3 days) C57BL/6 murine splenocytes were stained with Anti-CD86 antibody (STJA0000104) [Used at $0.2 \mu\text{g}/10^6$ cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.