

Anti-Adgre1 antibody [Cl:A3-1] {FITC} (STJA0000122)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	13733
Gene Symbol	Adgre1
UniProt ID	AGRE1_MOUSE

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Mouse abdominal macrophages elicited by starch broth were stained with Anti-F4/80 antibody (STJA0000122) [Used at 0.2 $\mu\text{g}/10^6$ cells dilution] (filled gray histogram). Unstained macrophages (blank black histogram) were used as control.