

## Anti-Adgre1 antibody [Cl:A3-1] {PerCP/Cyanine5.5} (STJA0000128)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rat                |
| Reactivity    | Mouse              |

### PRODUCT PROPERTIES


|                     |                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                       |
| Concentration       | 0.2 mg/mL                                                                                        |
| Conjugation         | PerCP/Cyanine5.5                                                                                 |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Isotype             | IgG2bk                                                                                           |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">13733</a>       |
| Gene Symbol | <a href="#">Adgre1</a>      |
| UniProt ID  | <a href="#">AGRE1_MOUSE</a> |

### ADDITIONAL INFORMATION

**Note** FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Mouse abdominal macrophages elicited by starch broth were stained with Anti-F4/80 antibody (STJA0000128) [Used at 0.2  $\mu$ g/ $10^6$  cells dilution] (filled gray histogram). Unstained macrophages (blank black histogram) were used as control.