

## Anti-CD62L/Sell antibody [Me114] {FITC} (STJA0000161)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rat                |
| Reactivity    | Mouse              |

### PRODUCT PROPERTIES


|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 5 Mu L/Test  |
| Conjugation         | FITC   |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Isotype             | IgG2ak   |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">20343</a>       |
| Gene Symbol | <a href="#">Sell</a>        |
| UniProt ID  | <a href="#">LYAM1_MOUSE</a> |

### ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



C57BL/6 murine splenocytes were stained with Anti-CD62L antibody (STJA0000161) (filled gray histogram). Unstained lymphocytes (empty black histogram) were used as control.