

## Anti-CD3 antibody [17A2] {FITC} (STJA0000181)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rat                |
| Reactivity    | Mouse              |

### PRODUCT PROPERTIES


|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 5 Mu L/Test  |
| Conjugation         | FITC   |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Isotype             | IgG2bk   |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|         |  |
|---------|--|
| Gene ID | <a href="#">12500</a><br><a href="#">12503</a><br><a href="#">12501</a><br><a href="#">Cd3d</a><br><a href="#">Cd247</a><br><a href="#">CD3D_MOUSE</a><br><a href="#">CD3Z_MOUSE</a><br><a href="#">CD3E_MOUSE</a> |
|---------|--|

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|



Mouse splenocytes were stained with Anti-CD3 antibody (STJA0000181) (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.