

Anti-CD3 antibody [17A2] {PE} (STJA0000182)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES


| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 5 Mu L/Test |
| Conjugation | PE |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2bk |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------|--|
| Gene ID | 12501 12503 12502 Cd3e Cd247 CD3E_MOUSE CD3Z_MOUSE CD3G_MOUSE |
|---------|--|

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine splenocytes were stained with Anti-CD3 antibody (STJA0000182) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.