

Anti-CD3 antibody [17A2] {PerCP} (STJA0000184)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PerCP
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID
[12501](#)
[12502](#)
[12500](#)
[Cd3e](#)
[Cd3g](#)
[CD3E_MOUSE](#)
[CD3G_MOUSE](#)
[CD3D_MOUSE](#)

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Mouse splenocytes were stained with Anti-Mouse CD3 Monoclonal Antibody (PerCP Conjugated) (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.